# Projekt: TCP/IP Datenserver

Projekttyp: Software

Dokumentationstyp: Kundeninformation

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#### DataServerTCP – Protocols : Overview (Version 1)

protocol type	protocol name	event for transmission to TCP/IP client
1	general information protocol (GIP)	The user has to set up the EEG machine for recording (number of electrodes, sampling frequency, patient's name, etc). After selecting a new session (record) the customer switches to the measuring mode (selects the tap "Measuring"). Every time the user switches to the measuring mode, the NEURO PRAX software sends out a general information protocol. Additionally, at the end of the recording an updated general information protocol will be sent.
2	marker name protocol (MNP)	This protocol contains the available marker names and values. It will be sent out to the TCP/IP client after switching from the tap "Patient" to the tap "Measuring" for a new session/recording.
3	impedance status protocol (ISP)	After preparation of electrode-to-skin-contact the user can check the impedances with the NEURO PRAX system. Every time the impedances has been checked, the NEURO PRAX system sends out the impedance status protocol with the updated impedance results.
4	data protocol (DP)	This protocol contains the recorded and/or processed data, the current sample index, a.s.o.
5	buffer overflow protocol (BOP)	Notification to the clients (buffer overflow at server).

**Similarities:**

character encoding:	Latin-1
field delimiter:	\$
protocol fields:	fixed size, each field contains one string
string alignment:	Left-aligned for text string filled with blanks Right-aligned for value strings starting with blanks
no. of fields:	not fixed, depends on the number of EEG channels and on the protocol type
EEG/EXG data format:	binary data, float precision (4 Bytes), little-endian byte ordering
protocol start:	neuroConn\$<protocolType>\$<protocolName>\$<protocolVersion>\$.....
protocol end:	end\$
channel / sample indexing:	starts with 0

## 3.2 General Information Protocol (Version 1)

parameter number	parameter name	parameter alias	size(bytes)	string	example
			incl. Delimiter \$	alignment	
1	neuroConn protocol	ProtocolStart	10	left	neuroConn
2	protocol type	ProtocolType	4	right	1
3	protocol name	ProtocolName	18	left	DataServerTCP-GIP
4	protocol version	ProtocolVersion	4	right	1
5	recording: file name	recFileName	19	left	20100329151422.EEG
6	recording: path name	recPathName	255	left	d:\neuroprax\datafiles
7	database: patient's name	DB_PatName	255	left	Mr. Public
8	database: patient's first name	DB_PatFirstName	255	left	John Q.
9	database: patient's birthday	DB_PatBirthday	11	left	2010-04-08
10	database: patient's identification number	DB_PatINR	255	left	101
11	name of the electrode setup	elecSetup	255	left	EEG-27-EP
12	sampling frequency	smpFreq	6	right	4000
13	name of selected algorithm	selAlgorithm	255	left	TMS
14	total number of channels	numChannels	5	right	4
15	number of EXG channels	numEXGchannels	5	right	4
16	channel names	ChannelNames	9 x numChannels	left	Fp1,Fp2,F7,F3
17	channel types	ChannelTypes	9 x numChannels	left	EEG,EEG,EEG,EEG
18	channel units	ChannelUnits	9 x numChannels	left	uV,uV,uV,uV
19	channel references	ChannelReference	9 x numChannels	left	GND,GND,GND,GND
20	protocol end	ProtocolEnd	4	left	end

\* birthday format:  
YYYY-MM-DD  
(similar to ISO 8601)

protocol size:           fixed:                           1616  
                          variable:               36 \* numChannels

### 3.3 Marker Name Protocol (Version 1)

parameter number	parameter name	parameter alias	size(bytes)	string	example
			incl. Delimiter \$	alignment	
1	neuroConn protocol	ProtocolStart	10	left	neuroConn
2	protocol type	ProtocolType	4	right	2
3	protocol name	ProtocolName	18	left	DataSetServerTCP-MNP
4	protocol version	ProtocolVersion	4	right	1
5	marker count	mC	4	right	10
6	marker index 1	markerIndex1	7	right	16384
7	marker name 1	markerName1	33	left	Start Record
...	...	...	...	...	...
...	marker index mC	markerIndex-mC	7	right	16384
...	marker name mC	markerName-mC	33	left	
...	protocol end	ProtocolEnd	4	left	end

protocol size:      fixed:                      44  
                          variable:                40 \* mC

### 3.4 Buffer Overflow Protocol (Version 1)

parameter number	parameter name	parameter alias	size(bytes)	string	example
			incl. Delimiter \$	alignment	
1	neuroConn protocol	ProtocolStart	10	left	neuroConn
2	protocol type	ProtocolType	4	right	5
3	protocol name	ProtocolName	18	left	DataSetServerTCP-BOP
4	protocol version	ProtocolVersion	4	right	1
5	protocol end	ProtocolEnd	4	left	end

protocol size:      fixed:                      40                      Bytes

### 3.5 Impedance Status Protocol (Version 1)

parameter number	parameter name	parameter alias	size(bytes) incl. Delimiter \$	string alignment	example
1	neuroConn protocol	ProtocolStart	10	left	neuroConn
2	protocol type	ProtocolType	4	right	3
3	protocol name	ProtocolName	18	left	DataSetTCP-ISP
4	protocol version	ProtocolVersion	4	right	1
5	number of EXG channels	numEXGchannels	5	right	4
6	channel name 0	chName0	9	left	Fp1
7	impedance status 0	impStatus0	3	right	0*
8	channel name 1	chName1	9	left	Fp2
9	impedance status 1	impStatus1	3	right	-1*
10	channel name 2	chName2	9	left	F7
11	impedance status 2	impStatus2	3	right	-2*
...					
...					
...					
(K-1)*2+6	channel name (K-1) **	chName(K-1)	9	left	O2
(K-1)*2+7	impedance status (K-1) **	impStatus(K-1)	3	right	-2*
(K-1)*2+8	protocol end	ProtocolEnd	4	left	end

\* impedance status

\*\*

-2 = poor contact, inacceptable (color: red) <as shown in the NEURO PRAX>

K denotes the number of EXG channels

-1 = worthy of improvement (color: yellow) <not shown in the NEURO PRAX>

0 = good contact (color: green) <as shown in the NEURO PRAX>

protocol size: fixed: 45  
variable: 12 \* numEXGchannels

### 3.6 Data Protocol (Version 1)

parameter number	parameter name	parameter alias	size(bytes) incl. Delimiter \$	string alignment	example
1	neuroConn protocol	ProtocolStart	10	left	neuroConn
2	protocol type	ProtocolType	4	right	4
3	protocol name	ProtocolName	18	left	DataServerTCP-DP
4	protocol version	ProtocolVersion	4	right	1
5	sample index	sampleIndex	12	right	2147483647
6	number of samples (=N)	numSmp	12	right	2147483647
7	number of EXG channels (=K)	numChannels	12	right	2147483647
8	data	data	<see below>	<binary>	
9	protocol end	ProtocolEnd	4	left	end

\*\*

The time stamp is measured in samples  
after start record.

#### data:

matrix data (=x) with K channels (K=numChannels) and N sample s(=numSmp)  
precision: float (4 Bytes)  
machine format: IEEE floating point with little-endian byte ordering

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- [1] Cheshire, S., Steinberg, D.H. Zero Configuration Networking. O'Reilly Media, Inc. 2006. ISBN 0-596-10100-7