


# NEURO PRAX<sup>®</sup> TMS/tES – Application Note

Part I: EEG & tDCS TCP/IP




# NEURO PRAX<sup>®</sup> TMS/tES – Application Note

## Part I: EEG & tDCS TCP/IP

neuroConn's full-band EEG system NEURO PRAX<sup>®</sup> TMS/tES measures EEG-Signals and Evoked Potential during:

- TMS (transcranial magnetic stimulation)
- tES (transcranial electric stimulation)




tDCS (transcranial direct current stimulation)  
tACS (transcranial alternating current stimulation)  
tRNS (transcranial random noise current stimulation)

The special NEURO PRAX<sup>®</sup> TMS/tES hard- and software offers:

- short recovery time 3-5 ms after TMS pulse
- high dynamic input range and resolution
- online artefact correction software for TMS and tACS artefacts
- full-band recording (0-1200 Hz) of EEG, EMG, ECG, EOG, ...
- online data interface to Brainsight<sup>®</sup> TMS navigation software
- export functions for analysis software (BESA, EMSE, ...)

# NEURO PRAX<sup>®</sup> TMS/tES – Application Note

Part I: EEG & tDCS TCP/IP



# NEURO PRAX<sup>®</sup> TMS/tES – Application Note

## Part I: EEG & tDCS TCP/IP

1. EEG electrodes from 32/ 64/ 128 electrodes and the special donut stimulation electrodes (anode and cathode)
2. The EEG electrodes are connected to the NEURO PRAX<sup>®</sup> TMS/tES input box of the amplifier.
3. The NEURO PRAX<sup>®</sup> TMS/tES amplifier is connected to the NEURO PRAX<sup>®</sup> TMS/tES panel PC.
4. The NEURO PRAX<sup>®</sup> TMS/tES panel PC is connected via a direct LAN cable or the same network to the external panel PC with Matlab.
5. In Matlab there is a special example to read data online from the NEURO PRAX<sup>®</sup> TMS/tES amplifier with a TCP/IP protocol.
6. In the Matlab computer you can program your special algorithm to detect e.g. slow wave sleep stages.
7. Once these stages are detected you can generate your particular stimulation parameters and output them with: sound card or DAC board.

# NEURO PRAX<sup>®</sup> TMS/tES – Application Note

## Part I: EEG & tDCS TCP/IP

8. This stimulation parameters can be entered into the DC-STIMULATOR PLUS with the REMOTE option. The DC-STIMULATOR PLUS will follow the particular given stimulation online.
9. The output current of the DC-STIMULATOR PLUS can be recorded with the NEURO PRAX<sup>®</sup> TMS/tES amplifier to monitor exactly how the stimulation is changing during the experiment. This is done using the SIGNAL OUT function and the SIGNAL OUT BOX.
10. In the case of tDCS and tACS, the NEURO PRAX<sup>®</sup> TMS/tES system can reject online the stimulation noise from the EEG data. Therefore, one can still monitor the actual EEG state during stimulation. During tDCS there is not special corrections but the amplifier will follow the offset and recover fast after the fade in. The high dynamic input range avoids saturation of the amplifier.
11. To be able to monitor each of the changes given a particular experiment or condition, the PC with matlab can send triggers via the trigger box to the NEURO PRAX<sup>®</sup> TMS/tES amplifier. These triggers will be used as markers in the EEG data for the different conditions used in the experiment.

# NEURO PRAX<sup>®</sup> TMS/tES – Application Note

Part I: EEG & tDCS TCP/IP

Contact



## neuroCare

Rindermarkt 7  
80331 München  
Germany  
Phone: +49 89 215 471 299 5

Albert-Einstein-Straße 3  
98693 Ilmenau  
Germany  
Phone: +49 3677 68 979 0

E-Mail: [info@neurocaregroup.com](mailto:info@neurocaregroup.com)  
Web: [www.neurocaregroup.com](http://www.neurocaregroup.com)



# NEURO PRAX<sup>®</sup> TMS/tES – Application Note

Part I: EEG & tDCS TCP/IP

Legal Notice

Document version: 1.0 (24/12/2016)

Copyright:

*All text, images, graphics and multimedia files and their layouts are subject to copyright and any other laws relating intellectual property. They may not be copied for commercial purposes or for transmission, modified or used in any other manner, in particular on other websites.*

*The content on this document is provided by neuroCare Group without assurance or warranty of any kind, whether explicitly or implicitly. All implicit guarantees concerning the merchantability, the suitability for certain purposes or the non-violation of laws and patents are also excluded. Even when we assume that the information that we provide is correct, it may still contain errors or inaccuracies.*

Trademarks:

*All trademarks mentioned on this website are the legally protected trade marks of the corresponding owners. They are mentioned for information purposes only. This applies in particular to model and brand names and corresponding logos and emblems of the relevant companies.*

Legal effect of this disclaimer:

*This disclaimer is considered to be part of the document which contained the link to this document. Should parts or individual formulations of this text not correspond to, no longer correspond to or incompletely correspond to the current legal situation, the validity of the content of the remaining parts of the document will be unaffected.*